

## DV1610 - ViewPoint Server Specifications

	Combined Application/Database Server (Minimum)	Combined Application/Database Server (Recommended)	Separate Application/Database Server (Minimum)	Separate Application/Database Server (Recommended)
<b>Operating Systems</b>	Windows Server 2008 R2	Windows Server 2012, 2012 R2	Windows Server 2008 R2	Windows Server 2012, 2012 R2
<b>Processor</b>	Intel 2.0 GHz Dual-Core or equivalent	Intel 2.6 GHz Quad-Core or equivalent	Intel 2.0 GHz Dual-Core or equivalent	Intel 2.6 GHz Quad-Core or equivalent
<b>RAM</b>	8 GB	16 GB	4 GB	8 GB
<b>HDD Disk Space</b>	60 GB	120 GB	30 GB	60 GB

**Note:** (1) With respect to all specifications in this document, hardware that is newer and higher-performance or a compatible software release that is newer is acceptable.

### Supported Web-Browsers:

- Microsoft Internet Explorer IE 10 or higher
- Microsoft Edge v20.X or higher
- Google Chrome v44.X or higher
- Safari v8 or higher
- Firefox v39.0 or higher
- Opera v32 or higher

### Additional Requirements:

- Microsoft .net 4.5
- Microsoft SQL Server (2008 R2, 2012, 2014 and 2016)
- ViewPoint Email Service enforces TLS encryption (TCP/IP port 587)

### Packages Installed by ViewPoint Installer:

- Java 7
- Ruby on Rails Framework
- TorqueBox
- Microsoft SQL Server 2014 Express